

Flex-T™ Tetramer and Cell Staining Protocol V.1

DOI

dx.doi.org/10.17504/protocols.io.e7bbhin



Kelsey Miller¹

¹BioLegend

BioLegend

Tech. support email: tech@biolegend.com



Sam Li

BioLegend



DOI: dx.doi.org/10.17504/protocols.io.e7bbhin

External link:

http://www.biolegend.com/media_assets/support_protocol/Protocol%20for%20fluorescent%20tetramer%20generation%20and%20cell%20staining%2006202016.pdf

Protocol Citation: Kelsey Miller . Flex-T™ Tetramer and Cell Staining Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e7bbhin>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: June 21, 2016

Last Modified: August 29, 2019

Protocol Integer ID: 3011

Abstract

Using UV-induced peptide exchange, MHC/peptide monomers can be generated with conditional Flex-T™ monomers that harbor peptides of interest in their binding grooves. These new MHC monomers are subsequently multimerized using streptavidin-fluorophore conjugates. The resulting Flex-T™ reagents can be used for staining antigen-specific T cells and flow cytometric analysis. In humans, the MHC molecules are called HLA (Human Leukocyte Antigen).

Guidelines

Materials

- Phosphate buffered saline pH 7.4, 10X concentrate (PBS, BioLegend Cat# 926201)
- Conditional Flex-T™ monomers
- 10 mM peptide solution of choice in 100% DMSO
- DMSO (e.g. Sigma-Aldrich Cat#D5879)
- 50 mM D-Biotin (e.g. Thermo Fisher, Cat#B20656)
- 10% (w/v) NaN₃ (e.g Sigma, Cat#S2002)
- Fluorophore-conjugated Streptavidin (BioLegend Cat# 405203, Cat# 405207, Cat#405225 or equivalent)
- Cell Staining Buffer (BioLegend Cat#420201 or equivalent)
- 96-well Polystyrene Microplate, U-shape (e.g. Falcon Cat#353077) or 5mL, 12 × 75mm tubes (e.g. Falcon Cat# 352008)
- Plate sealers (BioLegend Cat# 423601)
- 1.5 mL tubes (e.g. Eppendorf Cat# 022364111)

Equipment

- UV lamp, long-wave UV, 366 nm, 8 Watts (For example CAMAG cat# 022.9115, or Ultraviolet Crosslinker CL-1000)
- Incubator (37°C)
- Centrifuge capable of accommodating microtiter plates and tubes
- Single and multichannel pipettes capable of accurate delivery of variable volumes, and pipette tips

Precautions for use

- DMSO can be used to dissolve the peptides. However, do not exceed an end concentration of 10% (v/v) in the exchange reaction.
- Avoid repeated freeze-thawing.
- The Flex-T™/peptide solution needs to be kept on ice in the dark as much as possible. Do not work in front of a window.
- The use of short-wavelength (254 nm) or broad-band UV lamps is detrimental to MHC complexes.
- Centrifuge all vials before use (1 minute 3000xg at 4°C).

View protocol on Biolegend website for Representative Data



Peptide exchange

- 1 Bring all reagents to 0°C by putting them on ice.
- 2 Dilute 10mM stock solutions of peptides of choice to 400µM by mixing 5µl of peptide stock solution with 120µl PBS, and keep on ice.
- 3 Add 20µl diluted peptide (400µM) and 20µl conditional Flex-T™ monomer (200µg/mL) into 96-well U-bottom plate. Mix by pipetting up and down.
- 4 Seal the plate; centrifuge at 3300xg for 2 minutes at 4°C to collect the liquid down.
- 5 Remove the seal; put the plate on ice and illuminate with UV light for 30 minutes (the distance of the UV lamp to the samples should be 2-5 cm).
- 6 Seal the plate; incubate for 30 minutes at 37°C in the dark.
- 7 To evaluate the efficiency of the peptide exchange follow the Protocol for HLA class I ELISA to evaluate peptide exchange.

Generation of Tetramers

- 8 Transfer 30µl of peptide-exchanged monomer into a new plate, then add 3.3µl of conjugated streptavidin, mix by pipetting up-and-down. Incubate on ice in the dark for 30 minutes. This is enough for about 15 tests.

Note: BioLegend fluorophore-conjugated streptavidin products are recommended. For 30µl of exchanged Flex-T™ monomer we suggest to use 3.3µl of BioLegend PE-streptavidin (Cat#405203) or APC-streptavidin (Cat#405207). For BV421-streptavidin conjugate (Cat#405225) use 1.3µl. For optimal reaction with other fluorophore-conjugated streptavidin products ensure that the monomer:streptavidin conjugate has a 6:1 molar ratio.

- 9 During the incubation, prepare blocking solution by adding 1.6µl 50mM D-Biotin and 6µl 10% (w/v) Na₃N to 192.4µl PBS, mix by vortexing. After the incubation, add 2.4µl of blocking solution and pipette up and down to stop the reaction.
- 10 Seal the plate and incubate at 2 - 8°C overnight (or on ice for 30 minutes in the dark).



Note: We recommend Flex-T™ to be assembled with two different streptavidin conjugates in separate reactions. This allows for two-color staining with the same tetramer allele, ensuring the highest specificity.

Cell staining and flow cytometric analysis

- 11 Prepare cells of interest
- 12 Prior to perform staining, centrifuge the plate at 3300xg for 5 minutes at 4°C. Keep on ice in the dark.
- 13 Add 2×10^6 cells to a 96-well U-bottom plate or 12 × 75 mm tubes. Adjust volume to 200µl with CellStaining Buffer. Add 2µl per sample of Flex-T™ complex prepared in Steps 7 - 9, mix and incubate on ice in the dark for 30 minutes.
- 14 If co-staining with surface antibodies, prepare the antibody cocktail based on optimal staining concentration of each reagent. Incubate for 30 minutes on ice in the dark.
- 15 Wash the cells with Staining Buffer two times. Resuspend cells with Staining Buffer.
- 16 Acquire the samples with a flow cytometer and appropriate settings within 2 hours.

Note: A titration of the Flex-T™ is recommended for optimal performance.