

Mar 14, 2023

🌐 Calculating multiplicity of infection (MOI)

📖 [Arcadia](#)

DOI

dx.doi.org/10.17504/protocols.io.yxmvm2b2ng3p/v1

Januka Athukoralage¹, Adair Borges¹

¹Arcadia Science

Arcadia Science



Arcadia Science

Arcadia Science



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yxmvm2b2ng3p/v1

Protocol Citation: Januka Athukoralage, Adair Borges 2023. Calculating multiplicity of infection (MOI). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm2b2ng3p/v1>

Manuscript citation:

Athukoralage J, Borges A, Reiter T. (2023). Challenges of isolating bacteriophage mRNA for chemical analysis. <https://doi.org/10.57844/arcadia-j03a-mz33>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 09, 2023

Last Modified: January 29, 2024

Protocol Integer ID: 78460

Keywords: calculate, calculation, determine, estimate, establish, moi, multiplicity of infection, phage, phages, bacteriophage, bacteriophages, bacteria, infection, multiplicity, math, numbers, multiply, divide, colony, colonies, colony-forming units, CFU, CFUs, plaque, plaque-forming unit, plaque-forming units, unit, PFU, PFUs, culture, how to

Abstract

This is a simple protocol for the calculations that we routinely do to determine phage multiplicity of infection (MOI) in our experiments.

It's really important to know what MOI you are using. Ideally, you would perform these quantifications before your experiments such that you have a good ballpark understanding of what your multiplicity of infection is in your experiment. This also allows you to prepare phage dilutions to more precisely hit your target MOI.

We will also sometimes plate out bacteria and phage day-of using our experimental samples to give us an exact measurement of our MOI based on the culture that we used and to double-check that our phage dilution was correct.



Calculating colony-forming units per mL (CFU/mL)

- 1 Grow your bacteria to your target OD_{600} using the same culture conditions that you will use in your experiments.
- 2 Create a 10-fold dilution series of your bacterial culture. Use glass beads to seed an LB agar plate with 100 μ L of the 10^{-4} , 10^{-5} , 10^{-6} , and 10^{-7} dilutions of each subculture and incubate overnight. These dilutions typically yield a suitable number of colonies for accurate counting by eye. If they do not, you will need to adjust up or down depending on the exact details of what you are doing.

Count colonies the next day, and calculate how many colony-forming units (CFUs) are present in 1 mL of culture.

$$\text{CFU/mL} = (\# \text{ colonies} \times \text{dilution factor}) \div 0.1 \text{ mL}$$

Calculating plaque-forming units per mL (PFU/mL)

- 3 Infect 150 μ L bacterial culture with 10 μ L phage, respectively, from 10-fold dilution series of phage stocks. Aim for dilutions that will give you single countable plaques. Mix with with 3 mL of 0.5% (0.5 g/100 mL) agarose dissolved in LB supplemented with 1 mM $MgSO_4$, pour onto LB plates, and incubate at 37 °C overnight.

Count plaques the next day, and calculate how many plaque-forming units (PFUs) are present in 1 mL of phage stock.

$$\text{PFU/mL} = (\# \text{ plaques} \times \text{dilution factor}) \div 0.01 \text{ mL}$$

Calculating multiplicity of infection (MOI)

- 4 In our experiments, we generally aim for an MOI of ~ 10 to ensure that the majority of cells are infected with phage. However, you will need to determine the MOI appropriate for your experiments based on the biological phenomena you are studying. We also aim to have the phage volume be $\leq 10\%$ of the bacterial culture volume. So, we might use 500 μ L of bacterial culture and 50 μ L of phage stock. Again, you will need to determine which volumes are appropriate for your experiment. Here is how we calculate MOI using PFU/mL and CFU/mL.

$$\text{PFU/mL} \times \text{phage volume in mL} = \text{Total PFU}$$

$$\text{CFU/mL} \times \text{bacterial volume in mL} = \text{Total CFU}$$

Total PFU ÷ Total CFU = MOI