


Aug 15, 2022

DNA Extraction from Yeast

This protocol is a draft, published without a DOI.



Brian Teague¹

¹University of Wisconsin - Stout

Yeast ORFans CURE



Brian Teague

Trinity University

OPEN  ACCESS



Protocol Citation: Brian Teague 2022. DNA Extraction from Yeast. **Protocol exchange** <https://protocols.io/view/dna-extraction-from-yeast-cfagtibw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 15, 2022

Last Modified: August 15, 2022

Protocol Integer ID: 68648

Keywords: dna, extraction, saccharomyces, pcr, lioac, sds

Abstract

This is a "quick and dirty" way to get some genomic DNA out of yeast cells. Not pure enough for many things, but should be fine for a PCR template.

This protocol is adapted from

CITATION

Blount BA, Driessen MR, Ellis T (2016). GC Preps: Fast and Easy Extraction of Stable Yeast Genomic DNA.. Scientific reports.

LINK

<https://doi.org/10.1038/srep26863>

who in turn adapted it from

CITATION

Lõoke M, Kristjuhan K, Kristjuhan A (2011). Extraction of genomic DNA from yeasts for PCR-based applications.. BioTechniques.


LINK

<https://doi.org/10.2144/000113672>

Image Attribution

By gskx via Flickr. <https://www.flickr.com/photos/gskx/89462961>

Guidelines






The centrifugation steps all specify  21000 x g . If your microcentrifuge doesn't go this high, spin at the fastest speed available.

Materials

Equipment

- Dry bath at 72°C
- Dry bath at 42°C (optional)


Materials

- Sterile water
-  Lithium Acetate Dihydrate **Sigma Aldrich Catalog #L4158** solution, [M] 1 Molarity (M)
-  Sodium dodecyl sulfate **Sigma Aldrich Catalog #436143-25G** solution, [M] 10 Mass / % volume
-  Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500**
-  Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500** solution, [M] 70 % (v/v)
-  TE Buffer **Contributed by users**

Protocol materials

 Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500** In Materials, Materials and [2 steps](#)

 TE Buffer Materials, Step 10







 Lithium Acetate Dihydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L4158** Materials, Step 1

 Sodium dodecyl sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #436143-25G** Materials, Step 1

Safety warnings

- ! Both lithium acetate and SDS are irritants, particularly in the eyes. Wear appropriate PPE, including safety glasses, lab coats and gloves.

SDS is particularly gnarly if it's inhaled. If you're making a solution from powdered SDS, use a dust mask and/or weigh it out in a hood.



- 1 Make  100 μL of yeast lysis solution by mixing the following in a 1.7 ml microcentrifuge tube:
 -  70 μL H₂O
 -  20 μL  Lithium Acetate Dihydrate **Sigma Aldrich Catalog #L4158** solution, [M] 1 Molarity (M)
 -  10 μL  Sodium dodecyl sulfate **Sigma Aldrich Catalog #436143-25G** solution, [M] 10 Mass / % volume




- 2 Choose a yeast colony to analyze and circle it on the bottom of the petri dish. A large one is best.

- 3 Using a micropipette tip, scrape some of the colony off and resuspend it in the lysis solution. Vortex vigorously until the colony is mixed completely into the lysis solution.

Note

THIS IS NOT A CASE WHERE MORE IS BETTER. Your solution should be slightly cloudy. If it is quite "thick", then try again.




- 4 Incubate at  70 °C for  00:05:00 5m

- 5 Add  300 μL of  Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500** and vortex briefly. Centrifuge  21000 x g, 00:03:00 3m

Note

Make sure the centrifuge is balanced!

6 Using a P-1000 micropipettor, carefully aspirate the supernatant and discard as biological waste. **Do not disturb the pellet.**

7 Add  500 µL of [M] 70 % (v/v)  Ethanol (100%, Molecular Biology Grade) **Fisher Scientific Catalog #BP2818500** to the microcentrifuge tube. Centrifuge  21000 x g, 00:03:00 .

3m

Note


Make sure the centrifuge is balanced!



8 Using a P-1000 micropipettor, carefully aspirate the supernatant and discard as biological waste. **Do not disturb the pellet. Try to get as much of the ethanol off as you can.**


9 Let the pellet dry by leaving it in a  42 °C dry bath for  00:15:00 . Leave the cap open.

15m

Note

If you don't have a dry bath, just leave the tube (cap open) at room temperature. Extend the time to  00:30:00

10 Add  100 µL of  TE Buffer **Contributed by users** and vortex to resuspend the pellet.


11 Centrifuge  21000 x g, 00:00:30 to collect the cellular debris at the bottom. The DNA remains suspended in the supernatant.

30s

12 Label the tube and store at  -20 °C , or proceed directly to PCR.

13 If you need to use this sample again:

30s

- Thaw the sample completely.
- Vortex briefly to resuspend everything.
- Centrifuge  21000 x g, 00:00:30 to (re)collect the cellular debris at the bottom of the tube.

Citations

Blount BA, Driessen MR, Ellis T. GC Preps: Fast and Easy Extraction of Stable Yeast Genomic DNA.

<https://doi.org/10.1038/srep26863>

Lööke M, Kristjuhan K, Kristjuhan A. Extraction of genomic DNA from yeasts for PCR-based applications.

<https://doi.org/10.2144/000113672>