

Nov 21, 2019

## DNA/RNA Radiolabeling Protocol

DOI

[dx.doi.org/10.17504/protocols.io.8dshs6e](https://dx.doi.org/10.17504/protocols.io.8dshs6e)



Liz O'Brien<sup>1</sup>, Connor Tsuchida<sup>1</sup>

<sup>1</sup>University of California, Berkeley

The Center for Genome E...



Meredith Triplet

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8dshs6e](https://dx.doi.org/10.17504/protocols.io.8dshs6e)

**Protocol Citation:** Liz O'Brien, Connor Tsuchida 2019. DNA/RNA Radiolabeling Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8dshs6e>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 17, 2019

**Last Modified:** November 21, 2019

**Protocol Integer ID:** 28818

**Keywords:** Radiolabeling, DNA, RNA, CasX, TS, NTS

## Attachments



**Radiolabeling\_CasX D...**

141KB

## Guidelines

### CasX TS/NTS with non-hydrolysable spacers:

#### TS:

5'-CGCTAGCTACGT\***T**\***T**\***G**\***A**\***T**\***T**\***T**\***C**\***T**\***G**\***C**\***T**\***G**\***C**\***A**\***G**\***G**\***A**\*TGAAATCCCGTAATCGCGC-3'

MW: 15664.2 g/mol

Concentration: X  $\mu$ M

\*= phosphothioate, **bold letters** = PAM, *italic letters* = spacer

For 10 pmol of TS: X  $\mu$ l of stock

#### NTS:

5'-GCGCGATTACGGGAT**TT***CAT*\***C**\***C**\***T**\***G**\***C**\***A**\***G**\***C**\***A**\***G**\***A**\***A**\***A**\***A**\***T**\***C**\***A**\***A**\***A**\*CGTAGCTAGCG-3'

MW: 15749.3 g/mol

Concentration: X  $\mu$ M

\*= phosphothioate, **bold letters** = PAM, *italic letters* = spacer

For 10 pmol of NTS: X  $\mu$ l of substrate

### Labelling reaction setup:

#### \*TS:

XX  $\mu$ l DNA or RNA (10 pmoles)

2.5  $\mu$ l 10x PNK buffer

0.5  $\mu$ l PNK enzyme

1.5  $\mu$ l P32-gamma-ATP


XX mL dH<sub>2</sub>O (DEPC for labeling RNA) to 25  $\mu$ l


#### \*NTS:


XX  $\mu$ l DNA or RNA (10 pmoles)  
2.5  $\mu$ l 10x PNK buffer  
0.5  $\mu$ l PNK enzyme  
1.5  $\mu$ l P32-gamma-ATP  
XX mL dH<sub>2</sub>O (DEPC for labeling RNA) to 25  $\mu$ l


## Materials

### MATERIALS


 T4 Polynucleotide Kinase (3' phosphatase minus) - 200 units **New England Biolabs Catalog #M0236S**

 10X T4 PNK Reaction Buffer **New England Biolabs**

 ATP [ $\gamma$ -<sup>32</sup>P]- 3000Ci/mmol 10mCi/ml Lead 100  $\mu$ Ci (P32-gamma-ATP) **Perkin Elmer Catalog #NEG002A100UC**

 HiTrap Desalting columns with Sephadex G-25 resin **Ge Life Sciences Catalog #29048684**

## Safety warnings


 Please see SDS (Safety Data Sheet) for hazards and safety warnings.

1 Set up labeling reaction:



X $\mu$ l	DNA or RNA (10 pmoles)
2.5 $\mu$ l	10x PNK buffer
0.5 $\mu$ l	PNK enzyme
1.5 $\mu$ l	P32-gamma-ATP
	dH <sub>2</sub> O (DEPC for labeling RNA) to 25 $\mu$ l



Note

Mix the DNA, buffer, enzyme, and H<sub>2</sub>O at the bench, and then add the DNA/enzyme mixture to ATP-filled tubes in a radioactive use area.


2 Incubate at  37 °C for  00:30:00 .



3 Heat inactivate the PNK at  65 °C for  00:20:00 .

4 Prepare G25 columns (from GE, green box): vortex thoroughly, twist cap ¼ turn, snap off bottom, spin for  00:01:00 at  3000 rpm to get rid of liquid.






5 Add  50  $\mu$ L H<sub>2</sub>O to a labeled eppendorf tube, place G25 column in it.





6 Add  25  $\mu$ L H<sub>2</sub>O to each labeling reaction after heat inactivation is done.



7 Apply entire reaction (now 50  $\mu$ l total) to G25 column resin.

8 Spin for  00:02:00 at  3000 rpm . 

9 Since 50  $\mu$ l H<sub>2</sub>O were in bottom of tube and you add your 50  $\mu$ l reaction, you should end with up to 100  $\mu$ l of  100 nanomolar (nM) labeled DNA/RNA.

10 Measure  1  $\mu$ L of each reaction with the black rad counter on shelf to get cpm readings. 