


Feb 22, 2019

Glutaraldehyde fixation

DOI

dx.doi.org/10.17504/protocols.io.ycgfstw



Veronika Cencen¹

¹EPFL - EPF Lausanne



Veronika Cencen

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.ycgfstw

Protocol Citation: Veronika Cencen 2019. Glutaraldehyde fixation . [protocols.io https://dx.doi.org/10.17504/protocols.io.ycgfstw](https://dx.doi.org/10.17504/protocols.io.ycgfstw)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 20, 2019

Last Modified: February 22, 2019

Protocol Integer ID: 20584

Abstract

Fixation of adherent mammalian cells for AFM imaging

Guidelines

Prepare adherent cell culture on glass coverslips

Materials

MATERIALS

Glutaraldehyde, 25% solution **Bio Basic Inc. Catalog #GC3870.SIZE.500ml**

STEP MATERIALS

PBS

Protocol materials

Glutaraldehyde, 25% solution **Bio Basic Inc. Catalog #GC3870.SIZE.500ml** Materials

PBS Materials, Step 1

Safety warnings

! Glutaraldehyde requires a high degree of caution, even at low concentration! Carefully read SDS before use!

Classification according to Regulation (EC) No 1272/2008

Acute toxicity, Oral (Category 4), H302 Acute toxicity, Inhalation (Category 4), H332 Skin corrosion (Sub-category 1B), H314 Serious eye damage (Category 1), H318 Respiratory sensitisation (Category 1), H334 Skin sensitisation (Category 1), H317 Specific target organ toxicity - single exposure (Category 3), Respiratory system, H335 Short-term (acute) aquatic hazard (Category 1), H400 Long-term (chronic) aquatic hazard (Category 2), H411

Refer to SDS for details

Side effects include skin irritation.^[4] If exposed to large amounts, nausea, headache, and shortness of breath may occur.^[3] **Protective equipment** is recommended when used.^[3] Glutaraldehyde is effective against a range of **microorganisms** including **spores**.^{[3][6]}

As a strong sterilant, glutaraldehyde is **toxic** and a strong irritant.^[15] There is no strong evidence of **carcinogenic** activity.^[16] Some occupations that work with this chemical have an increased risk of some cancers.^[16]

Before start

Ensure proper disposal container is available in the area


Perform entire procedure in ventilated fumehood, disinfect workspace and tools with ethanol

Use gloves and eye protection (fumehood should be sufficient to prevent inhalation so a mask is not required)


- 1 Rinse cells 3 times with appropriate buffer, such as PBS or Ringer's

 PBS

- 2 Prepare concentration glutaraldehyde in KMgH buffer (or other buffer used for your application). Pipette just enough of the dilution on top of the rinsed cells to cover the area. Let sit for



Glutaraldehyde chemical formula



<http://publish.uwo.ca/~jkiernan/formglut.htm>



- 3 Rinse at least 3 times with buffer, properly disposing of the diluted glutaraldehyde.

Expected result

Glutaraldehyde is appropriate for whole and unroofed cell fixation, so up to 100% should be fixed with decent crosslinking. For thicker samples use an addition of paraformaldehyde

Safety information

Use sodium hydroxide for deactivation if necessary

http://msdssearch.dow.com/PublishedLiteratureDOWCOM/dh_0936/0901b80380936cd5.pdf?filepath=microbial/pdfs/noreg/253-03304.doc&fromPage=GetDoc