

Human T Cell Activation with anti-CD3 Antibodies (clone UCHT1, OKT3 or HIT3a) V.2

DOI

dx.doi.org/10.17504/protocols.io.hv7b69n



BioLegend, Inc.

BioLegend
Tech. support email: tech@biolegend.com



Kelsey Knight

BioLegend



DOI: dx.doi.org/10.17504/protocols.io.hv7b69n

External link: http://www.biolegend.com/media_assets/support_protocol/BioLegend_T_cell_Activation_Protocol_012715.pdf

Protocol Citation: BioLegend, Inc. . Human T Cell Activation with anti-CD3 Antibodies (clone UCHT1, OKT3 or HIT3a). [protocols.io](https://dx.doi.org/10.17504/protocols.io.hv7b69n)
<https://dx.doi.org/10.17504/protocols.io.hv7b69n>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: May 09, 2017

Last Modified: March 28, 2018

Protocol Integer ID: 5791

Guidelines

Materials:

- Sterile PBS
- Anti-human CD3 Antibody
 - Clone UCHT1 (LEAF™ format, Cat. No. 300413/300414/300432; Ultra-LEAF™ format, Cat. No. 300437/300438)
 - Clone OKT3 (LEAF™ format, Cat. No. 317303/317304/317315; Ultra-LEAF™ format, Cat. No. 317325/317326)
 - Clone HIT3a (LEAF™ format, Cat. No. 300313/300314; Ultra-LEAF™ format, Cat. No. 300331/300332)
- Cell culture medium (e.g., RPMI-1640 or IMDM supplemented with 10% FBS and 2mM L-glutamine)
- Sterile single-cell suspension of Ficoll-Hypaque-purified peripheral blood mononuclear cells, isolated T cells, or T cell subsets
- 96-well flat-bottom tissue culture plates with lids (e.g., Costar® Cat. No. 3596)

* Soluble forms of LEAF™ purified UCHT1 (1 µg/ml) or LEAF™ purified HIT3a (0.01 – 0.1 µg/ml) may be used to activate T cells from PBMC cell populations.

Materials

STEP MATERIALS

 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**

 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**

Protocol materials


 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**

In Materials, Materials, Step 2




- 1 Prepare a 10 µg/ml solution of anti-CD3 (clone UCHT1, OKT3, or HIT3a) in sterile PBS.

- 2 Dispense 50 µl of the antibody solution to each microwell of the 96-well assay plate. For the unstimulated control wells, add 50 µl of sterile PBS.

 96-well flat-bottom tissue culture plates with lids **Corning Catalog #3596**

- 3 Seal plate. Incubate at 37°C for 2 hours or 4°C overnight.

 02:00:00

- 4 Aseptically decant antibody solution from the microwell plate.

- 5 Wash plate microwells 3 times with sterile PBS (wash 1/3).

- 6 Wash plate microwells 3 times with sterile PBS (wash 2/3).

- 7 Wash plate microwells 3 times with sterile PBS (wash 3/3). Discard liquid.

- 8 Prepare single cell suspension of cells of interest in supplemented cell culture medium to 1-2 x 10⁶/ml.

- 9 Aliquot 200 µl cell suspension into plate microwells. Cover with lid. Incubate at 37°C in 5% CO₂ and 100% humidity for 3 days.