

May 14, 2018

## 🌐 LB (Luria-Bertani) agar

DOI

[dx.doi.org/10.17504/protocols.io.p43dqyn](https://dx.doi.org/10.17504/protocols.io.p43dqyn)

Katy M Monteith<sup>1</sup>

<sup>1</sup>University of Edinburgh

Vale Lab

Molecular Biology



Katy M Monteith

University of Edinburgh

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.p43dqyn](https://dx.doi.org/10.17504/protocols.io.p43dqyn)

**Protocol Citation:** Katy M Monteith 2018. LB (Luria-Bertani) agar. protocols.io <https://dx.doi.org/10.17504/protocols.io.p43dqyn>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 14, 2018

**Last Modified:** May 14, 2018

**Protocol Integer ID:** 12155

**Keywords:** Luria-Bertani, LB agar

## Abstract

This protocol describes how to make 1 L of Luria-Bertani agar (also known as LB agar) for the culturing of bacteria on plates.

## Guidelines

1 L of LB agar makes ~20 plates.

At 4 °C unused plates can be stored for around 2 weeks prior to use.

Once used, LB agar plates containing bacterial colonies can be stored in air-tight container/bag at 4 °C for a maximum of 2 weeks.

Prior to disposal, plates containing bacterial colonies must be wrapped in Parafilm and autoclaved.

## Materials


### MATERIALS

 Agar ash 2.0-4.5% **Sigma Aldrich Catalog #A7002**

 LB broth, Miller **Thermo Fisher Scientific Catalog #BP1426-500**

 LB agar, Miller **Thermo Fisher Scientific Catalog #BP1425-500**

## Safety warnings

 Read appropriate SDS for LB broth, Miller and/or LB agar, Miller before use.

Bottle can be extremely hot immediately following autoclave. Use thermal gloves for handling.

- 1 Add 25g LB broth, Miller and 15g agar (or 40g LB agar, Miller) per litre of water.
- 2 Mix well by inverting bottle several times until powder dissolved.
- 3 Sterilize by autoclaving for 20 min at 15 psi (1.05 kg/cm<sup>2</sup>) on liquid cycle.
- 4 Following autoclave (while media still liquid but cool enough to safely handle bottle) pour LB agar from bottle into plates within a laminar flow hood. Allow ~50 mL per plate.

**Note**

The LB agar bottle must only be opened and the plates poured in a laminar flow hood to prevent contamination.

- 5 Leave plates in laminar hood until agar sets (~30 mins).
- 6 When agar has set, replace lid, invert plates and store in an air tight bag at 4 °C.
- 7 Plates should be removed from 4 °C at least 1 h before use and placed in 30 °C incubator. To prevent condensation from lid dropping onto agar surface, wipe lids dry with Kimwipe and store plates inverted.

**Note**

Reminder: only remove lid from LB agar plate while in a laminar flow hood to prevent contamination.