

Aug 01, 2017

Lminex Milliplex Soluble Cytokine Receptor 13-plex

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.hvib64e

Troy Kemp¹, Ligia Pint¹

¹Human Papilloma Virus Immunology Laboratory, Leidos Biomedical Research, Inc., Frederick National Laboratory for Cancer Research



Alison Van Dyke

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.hvib64e

External link: <https://doi.org/10.1371/journal.pone.0182359>

Protocol Citation: Troy Kemp, Ligia Pint 2017. Lminex Milliplex Soluble Cytokine Receptor 13-plex. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.hvib64e>

Manuscript citation:

Dyke ALV, Kuhs KAL, Shiels MS, Koshiol J, Trabert B, Lofffield E, Purdue MP, Wentzensen N, Pfeiffer RM, Katki HA, Hildesheim A, Kemp TJ, Pinto LA, Chaturvedi AK, Safaeian M (2017) Associations between self-reported diabetes and 78 circulating markers of inflammation, immunity, and metabolism among adults in the United States. PLoS ONE 12(7): e0182359. doi: [10.1371/journal.pone.0182359](https://doi.org/10.1371/journal.pone.0182359)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 09, 2017

Last Modified: March 06, 2018

Protocol Integer ID: 5770

Abstract

Lminex Milliplex Soluble Cytokine Receptor 13-plex manufacturer's protocol

- 1 PREPARATION OF SAMPLES/REAGENTS FOR IMMUNOASSAY
- 2 Preparation of Serum/Plasma Thaw Time: Thaw the samples completely on ice, mix well by vortexing and centrifuge (10,000 rpm, 10 minutes, 4°C) prior to use in the assay to remove particulates. Serum/Plasma samples from normal subjects should be diluted 1:5 using the Serum Matrix provided in the kit as the sample diluent (20 µL sample mixed with 80 µL Serum Matrix). If samples require dilution beyond 1:5, continue to use the Serum Matrix as the sample diluent.
- 3 Preparation of Antibody-Immobilized Beads For individual vials of beads, sonicate each antibody-bead vial for 30 seconds then vortex for 1 minute. Add 60 µL from each antibody-bead vial to the Mixing Bottle and bring final volume to 3.0 mL with Bead Diluent. Vortex the mixed beads well. Unused portion may be stored at 2-8°C for up to one month. Example: When using 13 antibody-immobilized beads, add 60 µL from each of the 13 bead sets to the Mixing Bottle. Then add 2.22 mL Bead Diluent.
- 4 Preparation of Quality Controls Reconstitution Time: Before use, reconstitute Quality Control 1 and Quality Control 2 with 250 µL deionized water. Invert the vial several times to mix and vortex. Allow the vial to sit for 5-10 minutes and then transfer the controls to appropriately labeled polypropylene microfuge tubes. Unused portion may be stored at -20°C for up to one month.
- 5 Preparation of Wash Buffer Bring the 10X Wash Buffer to room temperature and mix to bring all salts into solution. Dilute 30 mL of 10X Wash Buffer with 270 mL deionized water. Store unused portion at 2-8°C for up to one month.
- 6 Preparation of Serum Matrix Reconstitution Time: Add 5.0 mL Assay Buffer to the bottle containing lyophilized Serum Matrix. Mix well. Allow at least 10 minutes for complete reconstitution. Leftover reconstituted Serum Matrix should be stored at -20°C for up to one month.
- 7 Preparation of Human Soluble Cytokine Receptor Standard Reconstitution Time: 1.) Prior to use, reconstitute the Human Soluble Cytokine Receptor Panel Standard with 250 µL deionized water to give STD7. Invert the vial several times to mix. Vortex the vial for 10 seconds. Allow the vial to sit for 5-10 minutes and then transfer the standard to an appropriately labeled polypropylene microfuge tube. This standard will be termed STD7; the unused portion may be stored at -20°C for up to one month.
- 8) Preparation of Working Standards. Label six polypropylene microfuge tubes STD6, STD5, STD4, STD3, STD2, and STD1. Add 150 µL of Assay Buffer to each of the six tubes. Prepare serial dilutions by adding 50 µL of STD7 reconstituted standard to the STD6 tube, mix well and transfer 50 µL of the STD6 standard to the STD5 tube, mix well and transfer 50 µL of the STD5 standard to the STD4 tube, mix well and transfer 50 µL of the STD4 standard to STD3 tube, mix well and transfer 50 µL of the STD3 standard to the STD2 tube and mix well, transfer 50 µL of

the STD2 standard to the STD1 tube and mix well. The 0 pg/mL standard (Background) will be Assay Buffer. Standard Volume of Deionized Water to Add (mL) Volume of Standard to Add Original (STD7) 2500 Standard Concentration (pg/ml) Volume of Assay Buffer to Add (mL) Volume of Standard to Add STD6 15050 mL of STD7 STD5 15050 mL of STD6 STD4 15050 mL of STD5 STD3 15050 mL of STD4 STD2 15050 mL of STD3 STD1 15050 mL of STD2 Reconstituted STD6 STD5 STD4 STD3 STD2 STD1 Standard (STD7) After dilution, each tube has the following concentrations for each analyte: Standard Tube # sIL-4R, sIL-6R, sRAGE, sTNF-R1, sTNF-R2 (pg/ml) sCD30, sGP130, sIL-1R1 (pg/ml) sEGFR, sIL-1RII, sVEGF-R1, sVEGF-R2, sVEGF-R3 (pg/ml) 1

- 9 2
- 10 41222
- 11 8
- 12 74883
- 13 3390.61,9534
- 14 31,562.57,81353,1256,25031,250612,50025,000125,000750,000100,000500,000 IMMUNOASSAY PROCEDURE Allow all reagents to warm to room temperature (20-25°C) before use in the assay. Run the standards, controls, and samples in duplicate. Set the filter plate on a plate holder at all times during reagent dispensing and incubation steps so that the bottom of the plate does not touch any surface.
- 15 Prewet the filter plate by pipetting 200 µL of Wash Buffer into each well of the Microtiter Filter Plate. Seal and mix on a plate shaker for 10 minutes at room temperature (20-25°C).
- 16 Remove Wash Buffer by vacuum. Blot excess Wash Buffer from the bottom of the plate with an absorbent pad or paper towels. 10) Add 25 µL of each Standard or Control into the appropriate wells. Add 25 µL Assay Buffer to the 0 pg/mL standard (Background).
- 17 Add 25 µL of Assay Buffer to the sample wells.



- 18 Add 25 μ L of the Serum Matrix solution to the background, appropriate standards, and control wells.
- 19 Add 25 μ L of Sample into the appropriate wells.
- 20 Vortex Mixing Bottle and add 25 μ L of the mixed Beads to each well. (Note: During addition of Beads, shake bead bottle intermittently to avoid settling.)
- 21 Seal the plate with a plate sealer, cover it with the lid. Wrap a rubber band around the plate holder, plate and lid and incubate with agitation on a plate shaker 2 hours at room temperature (20-25°C).
- 22 Gently remove fluid by vacuum.
- 23 Wash plate 2 times with 200 μ L/well of Wash Buffer, removing Wash Buffer by vacuum filtration between each wash. Blot excess Wash Buffer from the bottom of the plate with an absorbent pad or paper towels.
- 24 Add 25 μ L of Detection Antibodies into each well. (Note: Allow the Detection Antibodies to warm to room temperature prior to addition.)
- 25 Seal, cover with lid, and incubate with agitation on a plate shaker for 1 hour. DO NOT VACUUM AFTER INCUBATION.20) Add 25 μ L Streptavidin-Phycoerythrin to each well containing the 25 μ L of Detection Antibodies.
- 26 Seal, cover with lid and incubate with agitation on a plate shaker for 30 minutes at room temperature (20-25°C).
- 27 Gently remove all contents by vacuum.
- 28 Wash plate 2 times with 200 μ L/well Wash Buffer, removing Wash Buffer by vacuum filtration between each wash. Wipe any excess buffer on the bottom of the plate with a tissue.
- 29 Add 150 μ L of Sheath Fluid to all wells. Resuspend the beads on a plate shaker for 5 minutes.
- 30 Run plate on Luminex 100™ IS.



- 31 Save and analyze the data using Bio-Plex Manager software. EQUIPMENT SETTINGS
Events: 50, per bead region Sample Size: 100 μ L Gate Settings 4335 to 10,000 Time Out 60 seconds
QUALITY CONTROL The ranges for each analyte in Quality Control 1 and 2 are provided on the card insert or can be located at the MILLIPORE website www.millipore.com/techlibrary/index.do using the catalog number as the keyword.