


Feb 16, 2024

🌐 Measuring dopamine release in human-derived iPSCs using High-Performance Liquid Chromatography (HPLC) coupled to Electrochemical Detection (ECD)



📁 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.q26g7p2bkgwz/v1

Kaitlyn ML Cramb^{1,2,3}, Richard Wade-Martins^{1,2,3}, Stephanie J Cragg^{1,2,3}

¹Oxford Parkinson's Disease Center and Department of Physiology, Anatomy & Genetics, University of Oxford, Dorothy Crowfoot Hodgkin Building, South Parks Road, Oxford, OX1 3QU, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

³Kavli Institute for Nanoscience Discovery, University of Oxford, Dorothy Crowfoot Hodgkin Building, South Parks Road, Oxford, OX1 3QU, UK

Metabolomics Protocols & Workflows
Tech. support email: bbmisraccb@gmail.com



Cláudia C. Mendes

University of Oxford

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.q26g7p2bkgwz/v1

Protocol Citation: Kaitlyn ML Cramb, Richard Wade-Martins, Stephanie J Cragg 2024. Measuring dopamine release in human-derived iPSCs using High-Performance Liquid Chromatography (HPLC) coupled to Electrochemical Detection (ECD). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g7p2bkgwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 12, 2024

Last Modified: February 16, 2024

Protocol Integer ID: 95107

Keywords: Dopamine release, HPLC-ECD, human iPSCs

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020370

Abstract

This protocol allows for the detection of evoked or tonic dopamine release from human iPSC-derived dopamine neurons using High-Performance Liquid Chromatography (HPLC) coupled to electrochemical detection.

Guidelines

This assay can be performed on any mature iPSC-derived dopamine neurons (DANs), but success has been particularly found on iPSC-DANs that have been grown in culture for 50-90 days following a **modified Krik's protocol**. Optimal density found to be plated at 500 000 / well in a 24 well cell culture dish. Collected samples can be stored long-term at -80°C.

All steps are performed at Room Temperature (RT).

Materials

Reagents:

- **3,4-Dihydroxyphenylacetic Acid** (Sigma-Aldrich, CAT# 11569-25MG)
- **5-Hydroxyindole-3-acetic Acid** (Sigma Aldrich, CAT# H8876-5G)
- CaCl₂
- **Dopamine hydrochloride** (Sigma-Aldrich, CAT# H8502-5G)
- **Homovanilic acid** (Sigma-Aldrich, CAT# 69673 - 25MG)
- KCl
- **L-Norepinephrine hydrochloride** (Sigma-Aldrich, CAT# 74480-100MG)
- MgCl₂
- NaCl
- **Phosphate-buffered saline** (PBS), pH 7.4 (PBS) (Life Technologies, CAT# 10010056)
- Perchloric Acid (PCA)

Equipment:

- 4.6 x 150 mm Microsorb C18 reverse-phase column
- Decade II ECD
- Electrode (Antec Layden)

Making Ringers Buffer for evoked DA release:

Add each reagent to 500 mL dH₂O

- 1.2 mM CaCl₂
- 148 mM NaCl
- 0.85 mM MgCl₂
- 40 mM KCl
- pH 7.4, 300 osm

Making Ringers Buffer for tonic DA release:

Add each reagent to 500 mL dH₂O

- 1.2 mM CaCl₂
- 148 mM NaCl
- 0.85 mM MgCl₂
- 2.7 mM KCl
- pH 7.4, 300 osm

Making Mobile Phase Solution:

- 13% HPLC grade methanol
- 0.12 M Sodium phosphate monobasic dihydrate (NaH₂PO₄)
- 0.8 mM Ethylenediaminetetraacetic acid (EDTA)
- 0.5 mM 1-Octanesulphonic acid sodium salt (OSA)
- pH 4.6



Preparation

- 1 Add 1 μL of perchloric acid (PCA) to a labelled light-protected 1.5 mL Eppendorf tube for each sample.
- 2 Prepare Ringers Buffers for tonic or evoked release (see **Materials**). Store short-term at 4°C.

Dopamine (DA) Release Assay

- 3 Aspirate media from hiPSC-DANs plated as stated in guidelines.
- 4 Wash hiPSC-DANs one time with 100 μL phosphate buffered saline (PBS).
- 5 Aspirate PBS from hiPSC-DANs.
- 6 Add 100 μL of Ringer's Buffer for 5 minutes.
- 7 Collect Ringer's Buffer in pre-labelled tubes containing PCA.
- 8 Snap freeze on dry ice (optional).

Note

Samples can be stored at -80°C long-term or immediately processed for HPLC.

HPLC-ECD

- 9 If snap-frozen and stored at -80°C , thaw samples on ice.
- 10 Spin samples at 10 000 g for 10 minutes.



- 11 Run samples on HPLC column using a mobile phase running at 1 mL/min on a 4.6 x 150 mm Microsorb C18 reverse-phase column and Decade II ECD with a glassy carbon working electrode (Antec Layden) set at 0.7 V with respect to an Ag/AgCl reference electrode.
- 12 Calculate concentrate of dopamine released in sample compared to known standards.