


Sep 07, 2018 Version 2

Modified ZN Staining Protocol V.2

 Version 1 is forked from [mZN Staining Protocol](#)

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.tb2eiqe

Asar Khan¹, Sumaira Shams¹, Saima Khan¹, Muhammad Iftikhar Khan², Sardar Khan³, Abid Ali¹

¹Department of Zoology Abdul Wali Khan University Mardan, KPK Pakistan;

²Centre for Biotechnology and Microbiology University of Swat, Pakistan;

³Department of Environmental Science, University of Peshawar



Asar Khan

Department of Zoology Abdul Wali Khan University Mardan 2320...

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.tb2eiqe**

Protocol Citation: Asar Khan, Sumaira Shams, Saima Khan, Muhammad Iftikhar Khan, Sardar Khan, Abid Ali 2018. Modified ZN Staining Protocol . [protocols.io](#) **<https://dx.doi.org/10.17504/protocols.io.tb2eiqe>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 07, 2018

Last Modified: September 07, 2018

Protocol Integer ID: 15450

Abstract

The Modified Ziehl-Neelsen stain (mZN stain) is a type of differential bacteriological stain used to identify acid-fast organisms, mainly *Mycobacteria*. Acid fast organisms are those which are capable of retaining the primary stain when treated with an acid (*fast=holding capacity*). Members of the Actinomycetes, genus *Nocardia* (*N. brasiliensis* and *N. asteroides* are opportunistic pathogens) are partially acid-fast. Oocysts of coccidian parasites, such as *Cryptosporidium* and *Isospora*, are also acid-fast. Hence they can also be detected and identified through mZN staining procedure.

Materials

MATERIALS

⊗ Carbol-Fuchsin

⊗ Distilled Water

⊗ Methanol **Sigma Aldrich Catalog #M3641**

⊗ Disposable Latex Gloves, Medium, 100/Box **Bio Basic Inc. Catalog #GL002M.SIZE.1PK**

⊗ Methylene Blue **Gold Biotechnology Catalog #M-680**

⊗ Microscope slides

⊗ Compound Microscope

⊗ ethanol **BBI Biotech**

⊗ Acid Alcohol

STEP MATERIALS

⊗ Carbol-Fuchsin

⊗ Acid Alcohol

⊗ Methylene Blue **Gold Biotechnology Catalog #M-680**


Protocol materials


 Acid Alcohol In Materials, Materials, Step 6

 Methylene Blue **Gold Biotechnology Catalog #M-680** In Materials, Materials, Step 7

 Disposable Latex Gloves, Medium, 100/Box **Bio Basic Inc. Catalog #GL002M.SIZE.1PK** Materials

 Carbol-Fuchsin In Materials, Materials, Step 4

 Distilled Water Materials

 Microscope slides Materials

 ethanol **BBI Biotech** Materials


 Methanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3641** Materials

 Compound Microscope Materials

- 1 The stool sample was Spread evenly on the middle of the slide with constant rotational movement.




⌚ 00:10:00 (5 to 10 minutes) for rotational movement


 3 mg (Amount of stool sample)

2 The slides were then placed on a dryer with smeared surface upwards to air-dry them.

 60 °C

 00:10:00 minutes

3 The dried smear was fixed with absolute methanol.

 00:05:00 or (3-5 minutes)

4 Now, the Carbol-fuchsin solution was added to the slide to cover the whole smear.



⊗ Carbol-Fuchsin

⌚ 00:20:00 minutes

- 5 The slides were washed gently with tap water with the help of a dropper.

Safety information

Do not expose the slides to the high pressure of tap water directly, rather it will be better to use a dropper for washing the slides.

- 6 After washing the slide, decolorizer (Acid Alcohol) was added to the smear and the slide washed again with tap water.

⊗ Acid Alcohol

🧪 3 mL or 4-6 drops

- 7 Then the counter stain (Methylene Blue) was added and left for 5 minutes and then washed the slide with clean water.


⊗ Methylene Blue **Gold Biotechnology Catalog #M-680**

🕒 00:05:00 minutes wait for methylene blue

- 8 The back side of the slides were cleaned with a tissue paper and put in the draining rack to air-dry.

🕒 00:05:00 minutes, wait for slide to dry

- 9 The smear was examined with the help of a compound microscope with 40x and 100x (immersion oil lens) objective and scanned thoroughly for parasite identification.



Equipment

new equipment

NAME

Olympus

BRAND

CH20i

SKU

Biological microscope , Anti-fungus treated optics , Built to last- Superior build quality

SPECIFICATIONS

