

Jul 23, 2019 Version 1

🌐 NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB #E7546) V.1



DOI

dx.doi.org/10.17504/protocols.io.4nngvde

Menna Teffera¹, New England Biolabs¹

¹New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com



Isabel Gautreau

New England Biolabs

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4nngvde

External link: <https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale7546.pdf?rev=221a72d9f6624cfda4fc2fe760691f02>

Protocol Citation: Menna Teffera, New England Biolabs 2019. NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB #E7546). [protocols.io](https://dx.doi.org/10.17504/protocols.io.4nngvde) <https://dx.doi.org/10.17504/protocols.io.4nngvde>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 24, 2019

Last Modified: June 22, 2023

Protocol Integer ID: 25006

Abstract

The NEBNext Ultra II End Repair/dA-Tailing Module is optimized to convert 500 pg–1 µg of fragmented DNA to repaired DNA having 5′ phosphorylated, 3′ dA-tailed ends.

This module is part of the Ultra™ II workflow, and is optimized for use with the NEBNext™ Ultra II Ligation Module (NEB #E7595), for Illumina®-compatible library construction.

This module is also compatible with some Oxford Nanopore MinION™ workflows.

This module is designed for use with NEBNext Singleplex or Multiplex Oligos for Illumina (NEB #E7350, #E7335, #E7500, #E7600 or #E7535), NEBNext Ultra II Ligation Module (NEB #E7595), and NEBNext Ultra II Q5 Master Mix (NEB #M0544).

Kits that include reagents for every step in the Ultra II DNA library construction workflow are also available (NEBNext Ultra II DNA Library Prep Kit for Illumina (NEB [#E7645](#)) and NEBNext Ultra II DNA Library Prep with Sample Purification Beads (NEB [#E7103](#)).

Guidelines


Safe Stop Point: This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.

Caution: Signifies a step in the protocol that has two paths leading to the same point.

Color: A color listed before or after a reagent name indicates the cap color of the reagent to be added.

Materials

MATERIALS

 NEBNext Ultra II End Prep Reaction Buffer **New England Biolabs Catalog #E7647**

 NEBNext Ultra II End Prep Enzyme Mix **New England Biolabs Catalog #E7646**

Before start

Starting Material: 500 pg–1 µg fragmented DNA. We recommend that DNA be sheared in 1X TE. If the DNA volume post shearing is less than 50 µl, add 1X TE to a final volume of 50 µl. Alternatively, 10 mM Tris-HCl, pH 8.0 or 0.1X TE can be used.



NEBNext End Prep

- 1 Mix the following contents in a sterile nuclease-free tube:

Component	Volume
(green) NEBNext Ultra II End Prep Enzyme Mix	3 μ l
(green) NEBNext Ultra II End Prep Reaction Buffer	7 μ l
Fragmented DNA	50 μ l
Total Volume	60 μl

- 2 Set a 100 μ l or 200 μ l pipette to 50 μ l and then gently pipette the entire volume up and down at least 10 times to mix thoroughly. Perform a quick spin to collect all liquid from the sides of the tube.

Note

Note: It is important to mix well. The presence of a small amount of bubbles will not interfere with performance.

- 3 Place in a thermocycler, with the heated lid set to $\geq 75^{\circ}\text{C}$, and run the following program:

1h

🕒 00:30:00 at 🌡️ 20 °C

🕒 00:30:00 at 🌡️ 65 °C

Hold at 🌡️ 4 °C

Note

Safe Stop Point: If necessary, samples can be stored at -20°C ; however, a slight loss in yield ($\sim 20\%$) may be observed. We recommend continuing with adaptor ligation before stopping.

- 4 Proceed directly to NEBNext Ultra II Ligation Module NEB #E7595.