

Feb 01, 2023

PhageFISH detailed protocol

 [PLOS One](#)

 Peer-reviewed method

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.4r3l273wqg1y/v1

Line Jensen Ostenfeld¹, Saria Otani¹

¹DTU

PLOS ONE Lab Protocols
Tech. support email: plosone@plos.org



Saria Otani

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4r3l273wqg1y/v1

Protocol Citation: Line Jensen Ostenfeld, Saria Otani 2023. PhageFISH detailed protocol. [protocols.io](https://dx.doi.org/10.17504/protocols.io.4r3l273wqg1y/v1)
<https://dx.doi.org/10.17504/protocols.io.4r3l273wqg1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 27, 2023

Last Modified: February 01, 2023

Protocol Integer ID: 75967

Keywords: Staining and embedding, CARD amplification, Phage probe hybridisation, Antibody binding

Funders Acknowledgement:

NovoNordisk

Grant ID: NNF16OC0021856

Abstract

This protocol details about PhageFISh protocol.

Attachments



627-1301.docx

32KB

Guidelines

Controls to consider:

- Faecal sample with no target for the phage probe

Timeframe:

Day 1		3h20m
Prepare samples	30 minutes*	
Fix samples	1 hour + 10 min (1h incubation)	
<i>Prepare <u>permeabilisation buffer and HCl</u></i>	10-15 min	
<i>Prepare ice for <u>permeabilisation</u></i>	5 min	
Wash	5 min	❄
Permeabilise cells	1 hour + 10 min (1h incubation)	
Wash	10 min	
Inactivate peroxidases	15 min (10 min incubation)	
Wash	5 min	❄
Day 2		6h5m
Prepare probes	15 min	
Hybridisation of cyanine-labelled probes	3 hours + 10 min (3h incubation)	
Wash	20 min (15 min incubation)	❄
Pre-hybridisation of DIG-labelled probes	1 hour + 15 min (1h incubation)	
<i>Prepare probes</i>	20 min	
Hybridisation of DIG-labelled probes	1 hour + overnight (1h incubation)	🌙
Day 3		5h35m
Wash	2 hours + 15 min (30min + 1.5h incubation)	
<i>Prepare antibody washing and CARD buffers</i>	30 min	
Antibody binding	2 hours + 15 min (30min + 1.5h incubation)	
Wash	35 min	
CARD amplification	1 hour (45 min incubation)	
Wash	30 min	❄
Day 4		1-6h
Staining and sealing slides	1 hour*	❄
Microscopy	1-5 hours*	
Total:	17 hours (not incl. microscopy) (approx. 12 hours incubation time) 3-5 days	




*depending on number of samples

❄ Freezing and stopping possible after step

🌙 Overnight incubation after step

Materials

Necessary materials:

- Poly-L-lysine coated glass slides with writing area
- Pencil for writing (DO NOT use sharpie)
- Pipette tip lids for holding glass slides (one will fit four slides, collect one lid for each condition tested)
- Humidity chambers (one for each formamide concentration used simultaneously). Anaerobic growth chambers work well.
- Aluminium foil (to protect samples from light)
- Ice
- Fume hood
- Incubator set to  46 °C
- Incubator (or oven) set to  85 °C
- Water bath set to  48 °C
- Optimised and diluted Cy-labelled probes (see *Optimisation of formamide concentration*)
- Diluted phage probes (see *Buffers and Reagents*)
- All buffers (see *Buffers and Reagents*)
- Faecal samples of interest

Note

- If possible, samples should be submerged in plenty of buffer. Four slides can be submerged in 30-50ml in a pipette tip lid. For washing, very light agitation could be used (e.g. the shaking incubator set to 25rpm).
- For valuable solutions (like probe-solutions), only cover the sample area and handle with care. Use 500µl-1ml to cover sample area.
- All incubations are at room temperature unless specified.
- DO NOT allow samples to dry unless specified.
- When working with paraformaldehyde and formamide always work in the hood.
- After using humidity chambers, allow fumes to evaporate in fume hood overnight.



Before start


Prepare buffers (see **Preparation of Buffers for PhageFISH protocol**).



Fix faecal samples to glass slides

1 Mix a loopful faecal sample with  10-20 μL PBS (1X) and vortex thoroughly. 



2 Allow suspension to settle for  00:05:00 to avoid large debris. 

3 Take  10 μL of the supernatant and place on coated glass slide.

4 Smear the droplet thinly over the slide using a cover slip.


Note




Avoid smearing all the way to the edges.

5 Allow the sample to dry – this should not take more than  00:10:00 . 



Note

If not dry after 10 minutes, aspirate off excess liquid.


6 Work in fume hood. Overlay the slides with 1% paraformaldehyde (PFA). Ensure the whole sample area is covered (approx.  1 mL).

7 Incubate for  01:00:00 at  Room temperature in the fume hood. 

Note

This incubation should NOT exceed  01:00:00 ! 


8 Aspirate off excess PFA.

9 Wash in PBS for  00:01:00 .

1m

Note

If a lot of PFA remains on the sample, rinse twice in PBS.

FREEZING POINT – if necessary, samples can be rinsed in sterile water and 96% ethanol and air dried before freezing in closed box covered with aluminium foil at  -20 °C .





Permeabilise cells

10 Add lysozyme to permeabilisation buffer.



11 Overlay samples with permeabilisation buffer.

12 Incubate  On ice for  01:00:00 .

1h



13 Discard permeabilisation buffer.

14 Wash samples in PBS for  00:05:00 .

5m



15 Wash samples in sterile water for  00:01:00 .

1m



Inactivate peroxidases

16 Incubate samples in  0.01 Molarity (M) HCl for  00:10:00 .

10m



17 Wash samples in PBS for  00:05:00 .


5m



18 Wash samples in sterile water for  00:01:00 .

1m




19 Wash samples in 96% ethanol for  00:01:00 .

1m




20 Allow slides to dry on blotting paper or filter paper.

Note

FREEZING POINT – if necessary, samples can be frozen after drying. Store in closed container covered with aluminium foil at  -20 °C .

Cy-labelled probe hybridisation (16S rRNA probes)


21 Work in fume hood. Place a paper towel in the bottom of the hybridisation chamber and soak in formamide/milliQ solution corresponding to the hybridisation buffer concentration.



22 Overlay samples with hybridisation buffer-probe mix at  0.5 ng/μl of each probe and close humidity chamber.

23 Incubate at  46 °C for  03:00:00 .

3h



24 Prepare the washing buffer – heat to  48 °C .

25 Work in fume hood. Overlay the samples with washing buffer and incubate for  00:15:00 at  48 °C (in humidity chamber to avoid formamide fumes).

15m



26 Wash samples in sterile water.



27 Allow samples to dry.

Note

FREEZING POINT – if necessary, samples can be frozen after drying. Store in closed container covered with aluminium foil at $-20\text{ }^{\circ}\text{C}$.

Phage probe hybridisation

28 Work in fume hood. Place a paper towel in the bottom of the hybridisation chamber and soak in formamide/milliQ solution corresponding to the hybridisation buffer concentration.

29 Overlay samples with hybridisation buffer (no probes!) and close humidity chamber ($500\text{ }\mu\text{L}$ per slide).

30 Incubate for $01:00:00$ at $46\text{ }^{\circ}\text{C}$.

1h



31 Cover the samples with hybridisation buffer-probe mix at $10\text{ }\mu\text{g}/\mu\text{l}$ of each probe ($500\mu\text{l}$ per slide).

32 Place the dish back in the humidity chamber and incubate for $01:00:00$ at $85\text{ }^{\circ}\text{C}$.

1h



33 Immediately place the humidity chamber at hybridisation temperature **Overnight**.

1h



34 Wash slides.



34.1 Wash slides in gene washing buffer I for $00:01:00$. (1/3)

1m







34.2 Wash slides in gene washing buffer I for $00:01:00$. (2/3)




1m

- 34.3 Wash slides in gene washing buffer I for  00:01:00 . (3/3) 1m
- 34.4 Wash slides in gene washing buffer I for  00:30:00 at  42 °C . 30m
- 35 Wash slides. 
- 35.1 Wash slides in gene washing buffer II for  00:01:00 . (1/3) 1m
- 35.2 Wash slides in gene washing buffer II for  00:01:00 . (2/3) 1m
- 35.3 Wash slides in gene washing buffer II for  00:01:00 . (3/3) 1m
- 35.4 Wash slides in gene washing buffer II for  01:30:00 at  42 °C . 1h 30m
- 36 Wash slides in PBS for  00:01:00 . 1m






Antibody binding

- 37 Cover slides with antibody-blocking solution. Incubate for  00:30:00 . 30m

- 38 Discard antibody-blocking solution and cover with antibody binding solution. Incubate for  01:30:00 . 1h 30m

- 39 Wash slides. 
- 39.1 Wash slides in antibody washing solution for  00:01:00 . 1m

- 39.2 Wash slides in antibody washing solution for  00:10:00 . (1/3) 10m
- 39.3 Wash slides in antibody washing solution for  00:10:00 . (2/3) 10m
- 39.4 Wash slides in antibody washing solution for  00:10:00 . (3/3) 10m







CARD amplification

- 40 Mix  1 mL amplification buffer with  10 μ L H_2O_2 and  2 μ L Alexa tyramides (488). Vortex to mix. 
- 41 Cover slides with CARD buffer-tyramide mix (approx.  500 μ L per slide). Incubate at  37 $^{\circ}$ C for  00:45:00 . 45m 
- 42 Wash slides. 
- 42.1 Wash slides in PBS for  00:01:00 . 1m
- 42.2 Wash slides in PBS for  00:05:00 . 5m
- 42.3 Wash slides in PBS for  00:10:00 at  46 $^{\circ}$ C . 10m
- 42.4 Wash slides in PBS for  00:10:00 at  46 $^{\circ}$ C . 10m
- 43 Wash slides in sterile water for  00:01:00 . 1m 
- 44 Wash slides in 96% ethanol for  00:01:00 . 1m

Note

FREEZING POINT

Staining and embedding

- 45 Mix  1 mL SlowFade Gold antifade reagent with 1 5m/ml DAPI (final concentration  5 µg/mL), can be stored at  Room temperature). 
- 46 Place  10 µL solution in small droplets on the slides.
- 47 Place coverslip and press down gently to remove air pockets without disturbing the sample area.
- 48 Seal edges with clear nail polish.
- 49 Samples can now be stored at  -20 °C in covered container indefinitely.