

Feb 14, 2019 Version 3

## 🌐 Transfection of constructs in diplomemids to enhance the HR pathway, using RS-1 (3-(benzylamino) sulfonyl)-4-bromo-N-(4-bromophenyl) benzamide). V.3

DOI

[dx.doi.org/10.17504/protocols.io.x45fqy6](https://dx.doi.org/10.17504/protocols.io.x45fqy6)

Binnypreet Kaur<sup>1</sup>, Drahomíra Faktorová<sup>2</sup>, Julius Lukeš<sup>3</sup>

<sup>1</sup>[Institute of Parasitology Czech Academy of Sciences Faculty of Sciences University of South Bohemia Czech Republic ];

<sup>2</sup>[Institute of Parasitology Czech Academy of Sciences Faculty of Sciences University of South Bohemia Czech Republic];

<sup>3</sup>[Institute of Parasitology Czech Academy of Sciences Faculty of Sciences University of South Bohemia Czech Republic ] [Canadian Institute for Advanced Research Toronto Canada]

Protist Research to Opti...

Julius Lukes



**Binnypreet Kaur**

1 Institute of Parasitology, Biology Centre, Czech Academy o...

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.x45fqy6](https://dx.doi.org/10.17504/protocols.io.x45fqy6)

**Protocol Citation:** Binnypreet Kaur, Drahomíra Faktorová , Julius Lukeš 2019. Transfection of constructs in diplomemids to enhance the HR pathway, using RS-1 (3-(benzylamino) sulfonyl)-4-bromo-N-(4-bromophenyl) benzamide).. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x45fqy6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** February 14, 2019

**Last Modified:** February 14, 2019

Protocol Integer ID: 20349

## Abstract

RS-1 (RAD51-stimulatory compound 1) is a stimulator of homologous recombination (HR) protein RAD51



- 1 Step 1: Count the cells and plan to 1 nucleofection with  $5 \times 10^7$  cells for each construct.
  
- 2 Step 2: Harvest the cells by centrifugation at 1300xg for 5 min at room temperature in Swing Bucket Rotor.
  
- 3 Step 3: Resuspend the cell pellet in 100ul of AMAXA Human T- cell solution at 4C (from refrigerator combine 81.8ul of Human T-cell nucleofector solution + 18.2ul Supplement).
  
- 4 Step 4: Add 5-10ug of (PCR) (linearized DNA) along with 0.3 uM RS-1 enhancer into the cuvette (resuspend in 10ul of H<sub>2</sub>O).
  
- 5 Step 5: Put everything into the cuvette, close the cap and place in the electroporator, cuvette should only fit in one direction, but metal sides should face towards you.
  
- 6 Step 6: Press for the Program X-001 to electroporate and repeat the same as mentioned in transformation protocol.
  
- 7