

Feb 01, 2019

## Zymo Duet DNA/RNA MiniPrep Plus Extractions

DOI

[dx.doi.org/10.17504/protocols.io.umteu6n](https://dx.doi.org/10.17504/protocols.io.umteu6n)

Amy Lyden<sup>1</sup>

<sup>1</sup>CZ Biohub

Chan Zuckerberg Biohub



Amy Lyden

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.umteu6n](https://dx.doi.org/10.17504/protocols.io.umteu6n)

**Protocol Citation:** Amy Lyden 2019. Zymo Duet DNA/RNA MiniPrep Plus Extractions. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.umteu6n>

**Manuscript citation:**

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 13, 2018

**Last Modified:** February 01, 2019

**Protocol Integer ID:** 16787

## Abstract

Attached are detailed instructions for following the kit guidelines for the Zymo Duet DNA/RNA MiniPrep Plus.

 \_D7003\_ZR-Duet\_DNA-RNA\_MiniPr...


## Guidelines

Follow kit guidelines for best practices.

 \_D7003\_ZR-Duet\_DNA-RNA\_MiniPr...


## Materials

### MATERIALS

 ZR-Duet™ DNA/RNA MiniPrep Plus **Zymo Research Catalog #D7003**

Reagent in Kit	Quantity	Storage Temperature	Notes on Preparation
DNA/RNA Lysis Buffer	50 ml	Room Temp.	
DNA/RNA Prep Buffer	50 ml	Room Temp	
DNA/RNA Wash Buffer (concentrate)	24 ml	Room Temp	Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml DNA/RNA Wash Buffer concentrate before use
DNase/RNase-Free Water	10 ml	Room Temp	
DNase I (lyophilized)	1	Room Temp (lyophilized), -20 C (reconstituted)	Prior to use, reconstitute the lyophilized DNase I with 275 µl DNase/RNase-Free Water. Mix by inversion. Store aliquots at -20° C.
DNA Digestion Buffer	4 ml	Room Temp	
DNA/RNA Shield™ (2x concentrate)	25 ml	Room Temp	Dilute in DNase/RNase-free Water to obtain 1X for storage
PK Digestion Buffer	5 ml	Room Temp	
Proteinase K	20 mg	Room Temp (lyophilized), -20 C (reconstituted)	Prior to use, reconstitute the lyophilized Proteinase K with 1040 µl Proteinase K Storage Buffer. Vortex to dissolve. Store at -20° C
Spin-Away™ Filters	50	Room Temp	
Zymo-Spin™ IICG Columns	50	Room Temp	
Collection Tubes	3x50	Room Temp	
<b>Other Reagents and Equipment</b>			
Bead bashing tubes (ZR BashingBead™ Lysis Tubes or Omni 2.8mm ceramic bead bashing tubes or other)			
Qiagen Tissulyser II			
Filtered Pipette tips			
Ethanol			
Eppendorfs (1.5 mL)			
Microcentrifuge			
Eppendorf rack			

## Safety warnings





 Operate in an appropriate BSL level space for your sample type and suspected infectious agents.

## Before start





Ensure your samples have been stored at -80C, preferably in RNA/DNA shield. Avoid freeze thaws to ensure sample integrity. Identify your sample type and follow the 'sample preparation' section for that type.





## Sample Preparation for cultured isolates

1. Add  320  $\mu\text{L}$  1X DNA/RNA shield (2x concentrate diluted in nuclease-free water) to bead bashing tubes
2. Select colonies and add to DNA/RNA shield in bead bashing tube.
3. Tissue lyse on Qiagen Tissue lyser II  00:02:00 minutes at a frequency of 1/30 s.  
Alternatively, vortex at the highest frequency for 2 minutes.
4. Spin  00:01:00 at 14000RCF.
5. Add  320  $\mu\text{L}$  of DNA/RNA Lysis buffer. Proceed with Nucleic Acid Purification.





## Sample Preparation for liquids (ie, plasma, tracheal aspirate, CSF)

- 2 1. Add 1 part sample to 3 parts DNA/RNA lysis buffer (ie,  100  $\mu\text{L}$  sample to  300  $\mu\text{L}$  lysis buffer) in a bead bashing tube
2. Tissue lyse on Qiagen Tissue lyser II  00:01:00 minutes at a frequency of 1/30 s.  
Alternatively, vortex at the highest frequency for 1 minutes.
3. Spin  00:01:00 at 14000RCF. Proceed with Nucleic Acid Purification.


## Sample Preparation for samples in DNA/RNA shield (1x)

- 3 1. If not homogenized, tissue lyse on Qiagen Tissue lyser II in bead bashing tube for desired amount of time at a frequency of 1/30 s. Alternatively, vortex at the highest frequency.
2. Add one part DNA/RNA Lysis buffer (ie  200  $\mu\text{L}$  of sample in DNA/RNA shield to  200  $\mu\text{L}$  of lysis buffer). Proceed with Nucleic Acid Purification.




## Sample Preparation for whole blood




- 4 1. Add  200  $\mu\text{L}$  of 2X DNA/RNA shield to each  200  $\mu\text{L}$  of b2 of blood sample and mix thoroughly.
2. For every  400  $\mu\text{L}$  of reagent/blood mixture, add  8  $\mu\text{L}$  of reconstituted Proteinase K and mix thoroughly.



3. Incubate at room temperature for  00:30:00 minutes
4. Add an equal volume of isopropanol to the reagent/blood mixture and mix by vortexing. Proceed with Nucleic Acid Purification.

## Nucleic Acid Purification

- 5 Label a Spin-Away Filter (yellow) and the microcentrifuge tube for each sample.
- 6 Transfer up to  700  $\mu\text{L}$  of the sample into a Spin-Away™ Filter (yellow) in a Collection Tube and centrifuge  00:00:30 seconds at 16000 RCF. If you have >  700  $\mu\text{L}$ , you can reload the column and centrifuge again. Save the flow-through if you wish to proceed with RNA Purification. If not, discard flow through, and skip RNA Purification steps.



(**For whole blood ONLY**, after doing the above, discard flow-through. Add  200  $\mu\text{L}$  of DNA/RNA Lysis buffer directly to the filter matrix. Let stand for  00:05:00 minutes . Then centrifuge  00:00:30 seconds at 16000 RCF. **SAVE** this flow through as your RNA is in it)

- 7 Proceed using the yellow column with DNA Purification section, and using the flow-through with the RNA purification section. After that, you will continue with the Nucleic Acid Purification for both DNA and RNA in parallel.








## DNA Purification

- 8 DNA is bound to the yellow column. Transfer the Spin-Away™ Filter (yellow) into a new Collection Tube. Proceed with Nucleic Acid Purification.










## RNA Purification

- 9 Transfer flow-through into Eppendorf if necessary. Add an equal volume of ethanol (95-100%) to the flow-through and mix well.
- 10 Transfer up to  700  $\mu\text{L}$  of the sample with ethanol into a Zymo-Spin™ IICG Column (green) in a Collection Tube and centrifuge  00:00:30 seconds at 16000 RCF. Discard the




flowthrough. Repeat by reloading the column if you have more sample remaining.

- 11 Wash the column with  400  $\mu\text{L}$  DNA/RNA Wash Buffer and centrifuge  00:00:30 seconds at 16000 RCF. Discard the flowthrough.
- 12 Add 80  $\mu\text{L}$  DNase I Reaction Mix (  5  $\mu\text{L}$  reconstituted DNase +  75  $\mu\text{L}$  DNA Digestion Buffer) directly to the column matrix. (can make a master mix of the same proportions for multiple samples)
- 13 Incubate the column at room temperature (  20  $^{\circ}\text{C}$  -  30  $^{\circ}\text{C}$  ) for  00:15:00 minutes . Proceed with Nucleic Acid Purification.

## Nucleic Acid Purification

- 14 Add  400  $\mu\text{L}$  DNA/RNA Prep Buffer to your (yellow/green) column and centrifuge  00:00:30 seconds at 16000RCF. Discard the flow-through.
- 15 Add  700  $\mu\text{L}$  DNA/RNA Wash Buffer and centrifuge  00:00:30 seconds at 16000RCF. Discard the flow-through.
- 16 Add  400  $\mu\text{L}$  DNA/RNA Wash Buffer and centrifuge the column for  00:02:00 minutes at 16000RCF to ensure complete removal of the wash buffer. Carefully transfer the column into clean microcentrifuge/Eppendorf tubes.
- 17 Add  50  $\mu\text{L}$  DNase/RNase-Free Water directly to the column matrix, let stand for  00:05:00 minutes , and then centrifuge  00:00:30 seconds at 16000RCF to elute DNA and RNA from the respective column.



- 18 Take all  50  $\mu$ L of eluted DNA/RNA in microcentrifuge tube and reload onto column. Let stand for  00:05:00 minutes , and then centrifuge  00:01:00 minute at 16000RCF to elute more DNA and RNA from the respective column.
- 19 Aliquot RNA into single use aliquot tubes and freeze at -80. Freeze DNA at -20 or -80.