Sep 08, 2017



Real-time RT-PCR

Forked from Real-time PCR

PLOS One

DOI

dx.doi.org/10.17504/protocols.io.iencbde

Yuejun Wang¹

¹Department of Prosthodontics, Peking University School and Hospital of Stomatology

Yuejun Wang



DOI: dx.doi.org/10.17504/protocols.io.iencbde

External link: https://doi.org/10.1371/journal.pone.0184182

Protocol Citation: Yuejun Wang 2017. Real-time RT-PCR. protocols.io https://dx.doi.org/10.17504/protocols.io.iencbde

Manuscript citation:

Wang Y, Liu Y, Fan Z, Liu D, Wang F, Zhou Y (2017) IGFBP2 enhances adipogenic differentiation potentials of mesenchymal stem cells from Wharton's jelly of the umbilical cord *via* JNK and Akt signaling pathways. PLoS ONE 12(8): e0184182. doi: <u>10.1371/journal.pone.0184182</u>

License: This is an open access protocol distributed under the terms of the **<u>Creative Commons Attribution License</u>**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 12, 2017

Last Modified: January 16, 2018

Protocol Integer ID: 6318

Total RNA isolation

1 Total RNA was isolated from MSCs with TRIzol Reagent (Invitrogen, Carlsbad, CA, USA).

Reverse transcription

2 Reverse transcription reactions contained 2 µg RNA, random hexamers or oligo (dT), and reverse transcriptase, and were performed according to the manufacturer's protocol (Invitrogen, Carlsbad, CA, USA).

Real-time PCR

3 Real-time PCR was performed with QuantiTect SYBR Green PCR kit (Qiagen, Hilden, Germany) and an Icycler iQ Multi-color Real-time PCR Detection System. The reactions were run in duplicate with 1 µl of cDNA template in a 20 µl reaction volume with the program running at 95°C for 15 min, followed by 40 cycles of 94°C for 15 s, 60°C for 30 s, and 72°C for 30 s. The amplification specificity was confirmed by melting curve. The mRNA level of the target gene was acquired from the value of threshold cycle (Ct) as a relative level to that of GAPDH through the formula 2 -∆Ct (∆Ct = GAPDH Ct - gene of interest Ct). The efficiency of the primers was confirmed by sequencing the conventional PCR products before applying for real-time PCR.